

PATIENT TRANSFER SLIDE



EASYMOVE

(PATIENT MOVER FROM BED TO BED)



Applicable in large, medium and/or small hospitals, in ICU, Operation theatre, Emergency room and Patient ward, etc.



Usage:

According to the current requirement of modern nursing science, nursing involves a lot of physical work, which leads to backache or lumbar pain on the part of nurse. While moving the patient, particularly the serious patient or the patient who needs special care, between stretchers, patient beds, trolleys, inspection beds, X-ray fluoroscopy tables and operation tables, EASY-MOVE is easy to use.

NURSE WORK SIMPLENESS

PATIENT FEELING SAFE



Objectives of the new invention:

Easymove is simple to use, cost-effective and light-weighted .It can not only relieve the hardship of doctors and nurses when moving patients, but also relieve pains and adverse affects of patient when being moved.

Make nursing work become simple and easy. Relieve the pains of patients.





®

EASYMOVE

(PATIENT MOVER FROM BED TO BED)

2014 New Item Easymove is the best tool to move the patients among stretchers, operating tables, trolleys, sickbeds, CT tables etc. For moving and nursing, it can relieve the pain of the patients by providing the steady and safe movement, and also relieve the labor intensity of the nurses. **Any weight patient is available.**

Specifications:

Model: Fly-330

Red: 330R

Green: 330G

Blue: 330B

Size: 170cm×50cm×1.5cm

Weight: 3.32kg

Material: Nylon, PC board

Pressure and bending resistance, environmental protection, folding



Maintenance and sterilization

Easymove can be cleaned by wet clothes. Hard brush on the gray board is forbidden. The cover can be cleaned in water at 60°C and sterilized in any common antiseptic with no erosion.

No need to draw out easymove from trolley and X-Ray table. It can be penetrated by X-Ray without influence of radiation.

Attention

- ① Nurses should take charge of the way to operate Easymove expertly;
- ② Gap should not be left or over 15cm between the sickbed and the trolley;
- ③ The trolley must be fixed up in case of slipping while moving;
- ④ Do not pull or push the bed sheet strongly in case of the accident;
- ⑤ Do not use as a stretcher.



MERINO INTERNATIONAL

Ahead....in medical technology innovations

B-40, (G.F) Pal Mohan Plaza, 11/56, Main Desh Bandhu Gupta Road,
Karol Bagh, New Delhi - 110 005 (INDIA)
Tel +91 11 455 88229, 2352 0809 Tele Fax +91 11 2352 0809 Cell +91 999 999 9615
E-mail merinointernational@gmail.com Web. www.merinointernational.com
An ISO 9001:2008 Accredited Company



M4570910IN
ISO 9001:2008

