



Sterile Ultra Violet Cabinet

Features:-

- * Suitable for storage, Drying & Decontaminating Sensitive Products like Respiratory accessories, 30 hollow tubes, 16 breathing bags and instruments.
- * Fully Automatic two glass door for visual inspection with Microprocessor based Monitor.
- * Automatically turns off the heated air when door of cabinet is opened.
- * 4 Nos. Ultra violet lamp 30 watts emitting radiations in germicidal band of 253.7 for Sterilization.
- * Sterilization with UV & HEPA filter (0.3 um). Conforming to standard 209b for HEPA filter leak testing. With spun fiber Prefilter.
- * The sterilization cycle starts automatically, each time the doors are closed. Insulated return casing in wall for hot air re-circulation providing faster drying and saving electricity.
- * Micro processor based temperature controller with adjustable temperature for heat sensitive parts & varying load requirements.
- * Timer to select the sterilization time 1 to 999 minute with Alarms & Elapsed time counter for UV Sources. Time 1 to 6 hours.
- * External indicator lamp to show the functioning of cabinet.
- * Equipped with Twin Fan flow turbine ventilation system to create vertical and horizontal air circulation over Ultra Violet Sources & Heating Devices.
- * Remains on even when opened through a concealed sterilization process.
- * Aluminum Interior & Exterior. with Slide through tube holders & Nylon coated full size 9 trays which can be easily rearranged
- * Dimension 500 x 1200 x 2000 mm, with holder to hold 30 tubes at a time. Power consumption 30 watts each.
- * 220V 50/60 Hz. Conforming to International safety & performance Standards for Unit & filters.

UVSC Grab..ET Deluxe



Models	Capacity in Litres	Resettable shelves	
UVSC Grab..ET Deluxe 450	500	9	



MERINO INTERNATIONAL

Head... in medical technology innovations

B-40, (G.F.) Pal Mohan Plaza, 11/56, Main Desh Bandhu Gupta Road, Karol Bagh, New Delhi - 110 005 (INDIA)

Tel +91 11 455 88229, 2352 0809 Tele Fax +91 11 2352 0809 Cell +91 999 999 9615
E-mail merinointernational@gmail.com Web. www.merinointernational.com